

RAPID TESTS



Rapid, Simple, Reliable Results

Cardiac Markers • Drugs of Abuse • Faecal Proteins • Food Analysis • Gastro-Enteric Diseases • Parasitology
Pregnancy • Respiratory Diseases • Sexual Transmitted Diseases • Urinary Proteins



Mascia Brunelli

RAPID TESTS


CARD

Immunochromatographic, rapid, qualitative test on card format.

Visually reading result through use of coloured marker.

The addition of the specimen (serum, urine, pre-treated faeces) determines the appearance of the result in a few minutes.


Two methods are proposed: "competitive" and "non competitive"; the results are showed to side.



DIPSTICK

The tests on dip-stick are based on the analytical method of the sandwich test.

The reactive strip comes to contact with the specimen through immersion. The vertical flow determines the result in a few minutes.



SEMIQUANTITATIVE (A.F.I.A. Method)

Semi-quantitative determination of analytes through direct deposition of specimen as spot on the reactive strip. The adding of eluent allows to obtain a visible result in a few minutes.



PRODUCT	CAT. N°	TESTS	PRODUCT	CAT. N°	TESTS
CARDIAC MARKERS			MICROBIOLOGY: upper and lower respiratory tract diseases		
COMBI CARDIAC MARKERS	VS81470	3X10 tests	ADENOVIRUS	VC1009	25 tests
CRP-C REACTIVE PROTEIN	VT84000	25 tests	ADENOVIRUS + RSV	VC1019	25 tests
D-DIMER	VS81500	25 tests	INFLUENZA A+B	VC1012	25 tests
TROPONIN I	VS81465	25 tests	RSV	VC1006	25 tests
DRUGS OF ABUSE			STREPTOCOCCUS A	VQ81202/VQ 81210	20 tests/50 tests
ALCOHOL TEST	VU80600	20 tests	STREPTOCOCCUS B	VQ81305/VQ 81310	20 tests/50 tests
AMPHETAMINE	VU85501/VU85510	20 tests/75 tests	MICROBIOLOGY: sexual transmitted diseases		
ANTIDEPRESSIVE (TCA)	VU86201	20 tests	CHLAMYDIA TRACHOMATIS	VQ81406/VQ81411	20 tests/60 tests
BARBITURATES	VU85901/VU85910	20 tests/75 tests	GONORRHOEA	VQ85400	25 tests
BENZODIAZEPINE	VU85701	20 tests	SYPHILIS (Ab)	VQ83000	25 tests
BUPRENORPHIN	VU86101/VU86110	20 tests/75 tests	VAGINITIS	VQ81500	10 tests
CANNABINOIDS (MARIJUANA)	VU85201/VU85210	20 tests/75 tests	MICROBIOLOGY: food analysis		
COCAINE	VU85401/VU85410	20 tests/75 tests	E.COLI O157	VC1010	25 tests
EDDP	VU85120	20 tests	SALMONELLA	VQ84060	10 tests
PHENICLIDINE	VU85601	20 tests	L.MONOCYTOGENES	VQ84050	10 tests
KETAMINE	VU86300	20 tests	MICROBIOLOGY: others		
METHADONE	VU85801	20 tests	DENGUE (IgG/IgM)	VQ84000	10 tests
METHAMPHETAMINE	VU85101/VU85110	20 tests/75 tests	HEPATITIS A (Ag)	VR80200	25 tests
METHYLENEDIOSMINTHAPHETAMINE	VU85150	20 tests	HEPATITIS A (IgM)	VR82005	30 tests
OPIATES (Eroin, Morphine, Codein)	VU85301/VU85310	20 tests/75 tests	LACTOFERRIN	VT81600/VT81605	25 tests/25 tests
OXICODINE	VU85130	20 tests	LEGIONELLA PNEUMOPHILA	VQ84075/VQ84100	25 tests/25 tests
DRUGS MULTI TEST 6 parameters	VU85000/VU85001	10x6 tests/20x6 tests	MONONUCLEOSIS (IgM)	VQ82705	25 tests
DRUGS MULTI TEST 7 parameters	VU85007	10x7 tests	TUBERCULOSIS (Ab)	VQ81800	10 tests
DRUGS MULTI TEST 10 parameters	VU85010	10x10 tests	PARASITOLOGY		
PREGNANCY			CRYPTOSPORIDIUM	VC1005	25 tests
Beta-HCG MONOSTEP Test	VP80416/VP80417	20 tests/100 tests	CRYPTOSPORIDIUM-GIARDIA	VC1021	25 tests
Beta-HCG MONOSTEP Test	VP80418	50 tests	FILARIASIS (Ab)	VQ85300	25 tests
MICROBIOLOGY: gastro-enteric diseases			GIARDIA	VC1013	25 tests
ADENOVIRUS	VC1002/VC194020	25 tests/25 tests	LEISHMANIA	VQ85200	10 tests
ADENO+ROTAVIRUS	VC1004/VC194025	25 tests/25 tests	LEPTOSPIRA	VQ85100	10 tests
ASTROVIRUS	VC1020	25 tests	MALARIA	VQ81700/VQ81705	25 tests/25 tests
CAMPYLOBACTER SPECIES Ag	VC1007	25 tests	URINARY PROTEINS		
C.DIFFICILE TOXIN A+B	VC194055	25 tests	BENCE JONES (A.F.I.A.)	VT82900/VT82910	100 tests/20 tests
E.COLI O157	VC1011	25 tests	KAPPA E LAMBDA (A.F.I.A.)	VT82970/VT83900	2x10 tests/100 tests
HELICOBACTER PYLORI (Ag)	VC1110/VC1150	10/50 tests	MICROALBUMIN MONOSTEP	VT81802	25 tests
HELICOBACTER PYLORI (Ag)+CONTROLS	VC1150P	50 tests	FAECAL OCCULT BLOOD		
HELICOBACTER PYLORI (Ag)	VT82000/VT82001	10/50 tests	HAEMOGLOBIN Card	VT81500/VT81510	10 tests/50 tests
HELICOBACTER PYLORI (Ab) Card	VQ81650/VQ81660	20 tests/50 tests	HAEMOGLOBIN / TRANSFERRIN	VT81650	25 tests
L.MONOCYTOGENES	VQ84050	10 tests	GUAIACOL	VT81560	100 tests
ROTAVIRUS	VC1001/VC1940022	25 tests/25 tests			

COLOUR CODEX: **Black** non-competitive technique (card), **Red** competitive technique (card), **Blu** dipstick, **Green** semiquantitative.
All the tests have minimum shelf -life 15 months - Store at room temperature.

IVDs - CE



Mascia Brunelli

Mascia Brunelli S.p.A. - Viale Monza 272, 20128 Milan, Italy - Tel. 02 25 2091 - www.masciabrunelli.it